




## Introduction



FIGAERO Schematic  
[Aerodyne Research, Inc.]

### FIGAERO Features

- 2 separate inlets
  - gas and particle
- Coupled with CIMS
- Less fragmentation than AMS

### Thermal Decomposition of FIGAERO

- FIGAERO is ramped up to 200 °C
- Temperature as low as 200 °C can lead to thermal decomposition (Stark et al., 2017)
- Limited knowledge on potential impact of thermal decomposition

### Objective of the Study

- To fully characterize the mechanisms of thermal decomposition (e.g., decarboxylation, dehydration, etc.) in the FIGAERO using a suite of standard compounds with different functional groups
- To assess the impact of ramping rate on thermal decomposition

## Method

### FIGAERO Instrument

- Coupled with ToF-CIMS
- Reagent ion: Iodide
- Mode: Desorption mode

### Experiment


- Solvent: acetone
- [Standard] prepared: 0.05g/L
- Mass loaded: 100ng

### Standard Setting


- Ramping (15 min)
- Ramped up to ~194 °C from room temperature
- Average rate: 10 °C/min
- Soaking (15 min)
- At a max. temperature
- Cooling (10 min)

## Experimental Group


### Group 1: C<sub>4</sub> Compounds




### Group 3: Carboxylic Acids



### Group 4: Alcohols




### Group 2: Different # of C




### Compounds' Names

- Group 1**
- A Succinic acid  
B Malic acid  
C Oxaloacetic acid  
D Tartaric acid
- Group 2**
- E 5-oxoazaleic acid  
F 4-oxoheptanedioic acid  
G 1,3-acetonedicarboxylic acid
- Group 3**
- H 1,3,5-benzenetricarboxylic  
I 1,3,5-cyclohexanetricarboxylic  
J Tricarballic acid  
K 1,2,3,4-butanetetracarboxylic  
L Citric acid
- Group 4**
- M Erythritol  
N Xylitol  
O Mannitol  
P Sucrose

## Results I. Standard Chemical Compounds



### Group 1: C<sub>4</sub> Compounds




### Group 2: Different # of C


Compound E	Compound F	Compound G
• Negligible dehydration (2%)	• Negligible dehydration	• Decarboxylation
	• ~2%	• ~50% parent decomposed

Decrease in carbonyl proximity  
Increase in thermal decomposition

### Group 3: Carboxylic Acids




### Group 4: Alcohols



## Results II. Ramping Rate Experiment

### Ramping Rate Experiment with Citric acid



- Slowest ramping rate differs the most from the fastest and the normal ramping rates
- Slowest ramping rate leads to less fragmented decomposition products and least decomposition
- Fastest and normal ramping rates have miniscule difference

## Discussion & Conclusion

- **Decarboxylation vs. Dehydration**
  - Dehydration
    - Most compounds
  - Decarboxylation (Compounds C, G, and H)
    - Stable resonance structure with ≥ 2 carboxylic acids
    - Beta-carbonyl adjacent to two carboxylic acids
    - Presence of ring that can form resonance structure (i.e. benzene) with more than 2 carboxylic acids
- **Characterization of Thermal Decomposition**
  - Degree of Thermal Decomposition
    - Most direct measure: Parental fraction
    - FIGAERO-CIMS's intrinsic parameters as x-axis
    - Did not observe a clear trend between degree of thermal decomposition and Molecular weight, nO, nC, O/C, etc.
- **Potential Implication of Thermal Decomposition in Data Interpretation**
  - Most common compounds in the ambient air are alkanes, diacids, and terpenes (Kroll et al., 2011 & Khan et al., 2000)
    - Alkanes and diacids → No thermal decomposition
    - Terpenes → Minor decarboxylation (2%)
  - For more oxidized products with abundant carboxylic acids with alcohols or carbonyls, thermal decomposition must be factored in → Ongoing research

## References & Acknowledgement

- **References**
- Khan, F.I. and Ghoshal, A.K., 2000. Removal of volatile organic compounds from polluted air. *Journal of loss prevention in the process industries*, 13(6), pp.527-545.
  - Kroll, J. H., Donahue, N. M., Jimenez, J. L., Kessler, S. H., Canagaratna, M. R., Wilson, K. R., & Mysak, E. R. (2011). Carbon oxidation state as a metric for describing the chemistry of atmospheric organic aerosol. *Nature Chemistry*, 3(2), 133.
  - Stark, H., Yatavelli, R. L., Thompson, S. L., Kang, H., Krechmer, J. E., Kimmel, J. R., & Campuzano-Jost, P. (2017). Impact of thermal decomposition on thermal desorption instruments: advantage of thermogram analysis for quantifying volatility distributions of organic species. *Environmental science & technology*, 51(15), 8491-8500.
- **Acknowledgement**
- This research was funded by NSF CAREER AGS-1555034, NOAA NA18OAR4310112, and PURA from Georgia Tech
  - The FIGAERO-HR-ToF-CIMS was purchased through NSF Major Research Instrumentation (MRI) grant 1428738